



# NRIP3 Antibody, Biotin conjugated

Product Code	CSB-PA868283LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQ35
Immunogen	Recombinant Human Nuclear receptor-interacting protein 3 protein (1-241AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear receptor-interacting protein 3 (Sarcoma antigen NY-SAR-105), NRIP3, C11orf14
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NRIP3