



NSUN4 Antibody, Biotin conjugated

Product Code	CSB-PA850270LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96CB9
Immunogen	Recombinant Human 5-methylcytosine rRNA methyltransferase NSUN4 protein (26-284AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	5-methylcytosine rRNA methyltransferase NSUN4 (EC 2.1.1.-) (NOL1/NOP2/Sun domain family member 4), NSUN4
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NSUN4