



NUBPL Antibody, Biotin conjugated

Product Code	CSB-PA819450LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TB37
Immunogen	Recombinant Human Iron-sulfur protein NUBPL protein (39-319AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Iron-sulfur protein NUBPL (IND1 homolog) (Nucleotide-binding protein-like) (huInd1), NUBPL, C14orf127
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.