



NUDCD3 Antibody, HRP conjugated

Product Code	CSB-PA811605LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IVD9
Immunogen	Recombinant Human NudC domain-containing protein 3 protein (142-361AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NudC domain-containing protein 3, NUDCD3, KIAA1068
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NUDCD3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.