



# NUDT11 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA836216LD01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q96G61
<b>Immunogen</b>	Recombinant Human Diphosphoinositol polyphosphate phosphohydrolase 3-beta protein (1-164AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Diphosphoinositol polyphosphate phosphohydrolase 3-beta (DIPP-3-beta) (DIPP3-beta) (hDIPP3beta) (EC 3.6.1.52) (Diadenosine 5',5'''-P1,P6-hexaphosphate hydrolase 3-beta) (Diadenosine hexaphosphate hydrolase (AMP-forming)) (EC 3.6.1.60) (Nucleoside diphosphate-linked moiety X motif 11) (Nudix motif 11) (hAps1), NUDT11, APS1 DIPP3B
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Others
<b>Target Names</b>	NUDT11
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.