



NUDT11 Antibody, HRP conjugated

Product Code	CSB-PA836216LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96G61
Immunogen	Recombinant Human Diphosphoinositol polyphosphate phosphohydrolase 3-beta protein (1-164AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Diphosphoinositol polyphosphate phosphohydrolase 3-beta (DIPP-3-beta) (DIPP3-beta) (hDIPP3beta) (EC 3.6.1.52) (Diadenosine 5',5'''-P1,P6-hexaphosphate hydrolase 3-beta) (Diadenosine hexaphosphate hydrolase (AMP-forming)) (EC 3.6.1.60) (Nucleoside diphosphate-linked moiety X motif 11) (Nudix motif 11) (hAps1), NUDT11, APS1 DIPP3B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NUDT11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.