



NUP58 Antibody, Biotin conjugated

Product Code	CSB-PA887135LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BVL2
Immunogen	Recombinant Human Nucleoporin p58/p45 protein (211-480AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleoporin p58/p45 (58 kDa nucleoporin) (Nucleoporin-like protein 1), NUP58, KIAA0410 NUPL1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NUP58
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.