





NXNL2 Antibody

Product Code	CSB-PA726124LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VZ03
Immunogen	Recombinant Human Nucleoredoxin-like protein 2 protein (1-135AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleoredoxin-like protein 2 (Rod-derived cone viability factor 2) (RdCVF2), NXNL2, C9orf121
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NXNL2

Image



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA726124LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA726124LA01HU at dilution of 1:100