



OCIAD1 Antibody, Biotin conjugated

Product Code	CSB-PA882129LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NX40
Immunogen	Recombinant Human OCIA domain-containing protein 1 protein (1-245AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	OCIA domain-containing protein 1 (Ovarian cancer immunoreactive antigen domain containing 1) (Ovarian carcinoma immunoreactive antigen), OCIAD1, ASRIJ OCIA
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	OCIAD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.