



OCIAD1 Antibody, Biotin conjugated

Product Code	CSB-PA882129LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NX40
Immunogen	Recombinant Human OCIA domain-containing protein 1 protein (1-245AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	OCIA domain-containing protein 1 (Ovarian cancer immunoreactive antigen domain containing 1) (Ovarian carcinoma immunoreactive antigen), OCIAD1, ASRIJ OCIA
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	OCIAD1