



OSBPL2 Antibody, Biotin conjugated

Product Code	CSB-PA884442LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H1P3
Immunogen	Recombinant Human Oxysterol-binding protein-related protein 2 protein (251-480AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oxysterol-binding protein-related protein 2 (ORP-2) (OSBP-related protein 2), OSBPL2, KIAA0772 ORP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	OSBPL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.