



OSGIN2 Antibody, Biotin conjugated

Product Code	CSB-PA897460LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y236
Immunogen	Recombinant Human Oxidative stress-induced growth inhibitor 2 protein (1-250AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oxidative stress-induced growth inhibitor 2 (hT41), OSGIN2, C8orf1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	OSGIN2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.