



OTUD6B Antibody, Biotin conjugated

Product Code	CSB-PA847646LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N6M0
Immunogen	Recombinant Human Deubiquitinase OTUD6B protein (1-293AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	OTU domain-containing protein 6B (EC 3.4.19.12) (DUBA-5), OTUD6B, DUBA5
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	OTUD6B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.