



OXR1 Antibody, FITC conjugated

Product Code	CSB-PA854043LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N573
Immunogen	Recombinant Human Oxidation resistance protein 1 protein (401-600AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oxidation resistance protein 1, OXR1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	OXR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.