




PAGE1 Antibody

Product Code	CSB-PA13219A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75459
Immunogen	Recombinant Human P antigen family member 1 protein (1-146AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	P antigen family member 1 (PAGE-1) (AL5) (G antigen 9) (GAGE-9) (G antigen family B member 1) (Prostate-associated gene 1 protein), PAGE1, GAGE9 GAGEB1
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	PAGE1

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA13219A0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.