



PANK1 Antibody, Biotin conjugated

Product Code	CSB-PA017421LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TE04
Immunogen	Recombinant Human Pantothenate kinase 1 protein (462-582AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pantothenate kinase 1 (hPanK) (hPanK1) (EC 2.7.1.33) (Pantothenic acid kinase 1), PANK1, PANK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PANK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.