



PANK1 Antibody, Biotin conjugated

Product Code	CSB-PA017421LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TE04
Immunogen	Recombinant Human Pantothenate kinase 1 protein (462-582AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pantothenate kinase 1 (hPank) (hPank1) (EC 2.7.1.33) (Pantothenic acid kinase 1), PANK1, PANK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PANK1