



PBLD Antibody, Biotin conjugated

Product Code	CSB-PA017501LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30039
Immunogen	Recombinant Human Phenazine biosynthesis-like domain-containing protein (11-288AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phenazine biosynthesis-like domain-containing protein (EC 5.1.-.-) (MAWD-binding protein) (MAWDBP) (Unknown protein 32 from 2D-page of liver tissue), PBLD, MAWBP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PBLD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.