




# PBXIP1 Antibody

<b>Product Code</b>	CSB-PA856896HA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q96AQ6
<b>Immunogen</b>	Recombinant Human Pre-B-cell leukemia transcription factor-interacting protein 1 protein (361-600AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:500-1:2000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Pre-B-cell leukemia transcription factor-interacting protein 1 (Hematopoietic PBX-interacting protein), PBXIP1, HPIP
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	PBXIP1

## Image



Western blot

All lanes: PBXIP1 antibody at 4µg/ml + A549 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 81, 78, 73 kDa

Observed band size: 81 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.