



PBXIP1 Antibody, FITC conjugated

Product Code	CSB-PA856896HC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96AQ6
Immunogen	Recombinant Human Pre-B-cell leukemia transcription factor-interacting protein 1 protein (361-600AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pre-B-cell leukemia transcription factor-interacting protein 1 (Hematopoietic PBX-interacting protein), PBXIP1, HPIP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PBXIP1