



PCOLCE2 Antibody, HRP conjugated

Product Code	CSB-PA017632LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKZ9
Immunogen	Recombinant Human Procollagen C-endopeptidase enhancer 2 protein (24-415AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Procollagen C-endopeptidase enhancer 2 (Procollagen COOH-terminal proteinase enhancer 2) (PCPE-2) (Procollagen C-proteinase enhancer 2), PCOLCE2, PCPE2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PCOLCE2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.