



PDCD1 Antibody, Biotin conjugated

Product Code	CSB-PA619964LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15116
Immunogen	Recombinant Human Programmed cell death protein 1 protein (21-170AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Programmed cell death protein 1 (Protein PD-1) (hPD-1) (CD antigen CD279), PDCD1, PD1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PDCD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.