



PDE6D Antibody, Biotin conjugated

Product Code	CSB-PA526126LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43924
Immunogen	Recombinant Human Retinal rod rhodopsin-sensitive cGMP 3',5'-cyclic phosphodiesterase subunit delta protein (1-150AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Retinal rod rhodopsin-sensitive cGMP 3',5'-cyclic phosphodiesterase subunit delta (GMP-PDE delta) (Protein p17), PDE6D, PDED
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PDE6D