



PDE6D Antibody, FITC conjugated

Product Code	CSB-PA526126LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43924
Immunogen	Recombinant Human Retinal rod rhodopsin-sensitive cGMP 3',5'-cyclic phosphodiesterase subunit delta protein (1-150AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Retinal rod rhodopsin-sensitive cGMP 3',5'-cyclic phosphodiesterase subunit delta (GMP-PDE delta) (Protein p17), PDE6D, PDED
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PDE6D
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.