



PDHA1 Antibody, Biotin conjugated

Product Code	CSB-PA017715LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08559
Immunogen	Recombinant Human Pyruvate dehydrogenase E1 component subunit alpha, somatic form, mitochondrial protein (223-390AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pyruvate dehydrogenase E1 component subunit alpha, somatic form, mitochondrial (EC 1.2.4.1) (PDHE1-A type I), PDHA1, PHE1A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PDHA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.