



PDIA5 Antibody, Biotin conjugated

Product Code	CSB-PA614533LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14554
Immunogen	Recombinant Human Protein disulfide-isomerase A5 protein (22-262AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein disulfide-isomerase A5 (EC 5.3.4.1) (Protein disulfide isomerase-related protein), PDIA5, PDIR
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PDIA5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.