



PDIK1L Antibody, HRP conjugated

Product Code	CSB-PA854018LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N165
Immunogen	Recombinant Human Serine/threonine-protein kinase PDIK1L protein (1-341AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase PDIK1L (EC 2.7.11.1) (PDLIM1-interacting kinase 1-like), PDIK1L, CLIK1L
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PDIK1L