



PECR Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA017769LA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9BY49 |
| Immunogen | Recombinant Human Peroxisomal trans-2-enoyl-CoA reductase protein (1-303AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Peroxisomal trans-2-enoyl-CoA reductase (TERP) (EC 1.3.1.38) (2,4-dienoyl-CoA reductase-related protein) (DCR-RP) (HPDHase) (Short chain dehydrogenase/reductase family 29C member 1) (pVI-ARL), PECCR, SDR29C1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cardiovascular |
| Target Names | PECCR |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |