



PECR Antibody, FITC conjugated

Product Code	CSB-PA017769LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BY49
Immunogen	Recombinant Human Peroxisomal trans-2-enoyl-CoA reductase protein (1-303AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peroxisomal trans-2-enoyl-CoA reductase (TERP) (EC 1.3.1.38) (2,4-dienoyl-CoA reductase-related protein) (DCR-RP) (HPDHase) (Short chain dehydrogenase/reductase family 29C member 1) (pVI-ARL), PECR, SDR29C1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PECR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.