



# PEX11B Antibody, HRP conjugated

|                     |                                                                                                                                      |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| Product Code        | CSB-PA017796LB01HU                                                                                                                   |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                        |
| Uniprot No.         | O96011                                                                                                                               |
| Immunogen           | Recombinant Human Peroxisomal membrane protein 11B protein (9-228AA)                                                                 |
| Raised In           | Rabbit                                                                                                                               |
| Species Reactivity  | Human                                                                                                                                |
| Tested Applications | ELISA                                                                                                                                |
| Form                | Liquid                                                                                                                               |
| Conjugate           | HRP                                                                                                                                  |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4                                                     |
| Purification Method | >95%, Protein G purified                                                                                                             |
| Isotype             | IgG                                                                                                                                  |
| Clonality           | Polyclonal                                                                                                                           |
| Alias               | Peroxisomal membrane protein 11B (Peroxin-11B) (Peroxisomal biogenesis factor 11B) (Protein PEX11 homolog beta) (PEX11-beta), PEX11B |
| Immunogen Species   | Homo sapiens (Human)                                                                                                                 |
| Research Area       | Tags & Cell Markers                                                                                                                  |
| Target Names        | PEX11B                                                                                                                               |