



PGF Antibody

Product Code	CSB-PA06959A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49763
Immunogen	Recombinant Human Placenta growth factor protein (19-170AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Placenta growth factor (PIGF), PGF, PGFL PLGF
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PGF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.