



PGK2 Antibody, HRP conjugated

Product Code	CSB-PA017859LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07205
Immunogen	Recombinant Human Phosphoglycerate kinase 2 protein (101-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphoglycerate kinase 2 (EC 2.7.2.3) (Phosphoglycerate kinase, testis specific), PGK2, PGKB
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PGK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.