



PHF10 Antibody, FITC conjugated

Product Code	CSB-PA017898LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WUB8
Immunogen	Recombinant Human PHD finger protein 10 protein (151-360AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PHD finger protein 10 (BRG1-associated factor 45a) (BAF45a) (XAP135), PHF10, BAF45A
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PHF10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.