





PI3 Antibody

Product Code	CSB-PA017952HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19957
Immunogen	Recombinant Human Elafin protein (61-117AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Elafin (Elastase-specific inhibitor) (ESI) (Peptidase inhibitor 3) (PI-3) (Protease inhibitor WAP3) (Skin-derived antileukoproteinase) (SKALP) (WAP four-disulfide core domain protein 14), PI3, WAP3 WFDC14
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	PI3

Image



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA017952HA01HU at dilution of 1:100



Immunofluorescent analysis of PC-3 cells using CSB-PA017952HA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)