



PI3 Antibody, Biotin conjugated

Product Code	CSB-PA017952HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19957
Immunogen	Recombinant Human Elafin protein (61-117AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Elafin (Elastase-specific inhibitor) (ESI) (Peptidase inhibitor 3) (PI-3) (Protease inhibitor WAP3) (Skin-derived antileukoproteinase) (SKALP) (WAP four-disulfide core domain protein 14), PI3, WAP3 WFDC14
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	PI3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.