



PIGO Antibody, Biotin conjugated

Product Code	CSB-PA855068LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TEQ8
Immunogen	Recombinant Human GPI ethanolamine phosphate transferase 3 protein (219-420AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GPI ethanolamine phosphate transferase 3 (EC 2.-.-.-) (Phosphatidylinositol-glycan biosynthesis class O protein) (PIG-O), PIGO
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PIGO
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.