



PIH1D1 Antibody, FITC conjugated

Product Code	CSB-PA889130LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NWS0
Immunogen	Recombinant Human PIH1 domain-containing protein 1 protein (1-290AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PIH1 domain-containing protein 1 (Nucleolar protein 17 homolog), PIH1D1, NOP17
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PIH1D1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.