



PIN1 Antibody, Biotin conjugated

Product Code	CSB-PA04619D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13526
Immunogen	Recombinant Human Peptidyl-prolyl cis-trans isomerase NIMA-interacting 1 protein (1-163AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidyl-prolyl cis-trans isomerase NIMA-interacting 1 (EC 5.2.1.8) (Peptidyl-prolyl cis-trans isomerase Pin1) (PPIase Pin1) (Rotamase Pin1), PIN1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PIN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.