



PITPNC1 Antibody, Biotin conjugated

Product Code	CSB-PA892330LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKF7
Immunogen	Recombinant Human Cytoplasmic phosphatidylinositol transfer protein 1 protein (1-268AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytoplasmic phosphatidylinositol transfer protein 1 (Mammalian rdgB homolog beta) (M-rdgB beta) (MrdgBbeta) (Retinal degeneration B homolog beta) (RdgBbeta), PITPNC1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PITPNC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.