



# PLA2G16 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA018089LD01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P53816
<b>Immunogen</b>	Recombinant Human HRAS-like suppressor 3 protein (1-132AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	HRAS-like suppressor 3 (HRSL3) (EC 3.1.1.32) (EC 3.1.1.4) (Adipose-specific phospholipase A2) (AdPLA) (Group XVI phospholipase A1/A2) (H-rev 107 protein homolog) (H-REV107) (HREV107-1) (HRAS-like suppressor 1) (HREV107-3) (Renal carcinoma antigen NY-REN-65), PLA2G16, HRASLS3 HREV107
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cancer
<b>Target Names</b>	PLA2G16
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.