



PLCL2 Antibody, Biotin conjugated

Product Code	CSB-PA892366LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UPR0
Immunogen	Recombinant Human Inactive phospholipase C-like protein 2 protein (2-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inactive phospholipase C-like protein 2 (PLC-L(2)) (PLC-L2) (Phospholipase C-L2) (Phospholipase C-epsilon-2) (PLC-epsilon-2), PLCL2, KIAA1092 PLCE2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PLCL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.