



PLEKHA8 Antibody, HRP conjugated

Product Code	CSB-PA856957LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96JA3
Immunogen	Recombinant Human Pleckstrin homology domain-containing family A member 8 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pleckstrin homology domain-containing family A member 8 (PH domain-containing family A member 8) (Phosphatidylinositol-four-phosphate adapter protein 2) (FAPP-2) (Phosphoinositol 4-phosphate adapter protein 2) (hFAPP2) (Serologically defined breast cancer antigen NY-BR-86), PLEKHA8, FAPP2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PLEKHA8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.