



PLEKHS1 Antibody, Biotin conjugated

Product Code	CSB-PA716341LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5SXH7
Immunogen	Recombinant Human Pleckstrin homology domain-containing family S member 1 protein (181-462AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pleckstrin homology domain-containing family S member 1 (PH domain-containing family S member 1) (Epididymis luminal protein 185) (hEL185), PLEKHS1, C10orf81
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PLEKHS1