



# PLEKHS1 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA716341LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q5SXH7
<b>Immunogen</b>	Recombinant Human Pleckstrin homology domain-containing family S member 1 protein (181-462AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Pleckstrin homology domain-containing family S member 1 (PH domain-containing family S member 1) (Epididymis luminal protein 185) (hEL185), PLEKHS1, C10orf81
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Others
<b>Target Names</b>	PLEKHS1