



PLEKHS1 Antibody, HRP conjugated

Product Code	CSB-PA716341LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5SXH7
Immunogen	Recombinant Human Pleckstrin homology domain-containing family S member 1 protein (181-462AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pleckstrin homology domain-containing family S member 1 (PH domain-containing family S member 1) (Epididymis luminal protein 185) (hEL185), PLEKHS1, C10orf81
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PLEKHS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.