



POLR2E Antibody, HRP conjugated

Product Code	CSB-PA018332LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19388
Immunogen	Recombinant Human DNA-directed RNA polymerases I, II, and III subunit RPABC1 protein (1-210AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA-directed RNA polymerases I, II, and III subunit RPABC1 (RNA polymerases I, II, and III subunit ABC1) (DNA-directed RNA polymerase II 23 kDa polypeptide) (DNA-directed RNA polymerase II subunit E) (RPB5 homolog) (XAP4), POLR2E
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	POLR2E
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.