




POMGNT1 Antibody

Product Code	CSB-PA855513LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WZA1
Immunogen	Recombinant Human Protein O-linked-mannose beta-1,2-N-acetylglucosaminyltransferase 1 protein (59-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein O-linked-mannose beta-1,2-N-acetylglucosaminyltransferase 1 (POMGnT1) (EC 2.4.1.-) (UDP-GlcNAc:alpha-D-mannoside beta-1,2-N-acetylglucosaminyltransferase I.2) (GnT I.2), POMGNT1, MGAT1.2
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	POMGNT1

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA855513LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.