




PORCN Antibody

Product Code	CSB-PA887958LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H237
Immunogen	Recombinant Human Protein-serine O-palmitoleoyltransferase porcupine protein (269-332AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein-serine O-palmitoleoyltransferase porcupine (EC 2.3.1.-) (Protein MG61), PORCN, MG61 PORC PPN
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	PORCN

Image



Western blot


All lanes: PORCN antibody at 3µg/ml + Mouse thymus tissue

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 53, 52, 43 kDa

Observed band size: 53 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA887958LA01HU at dilution of 1:100



Immunofluorescence staining of A549 cells with CSB-PA887958LA01HU at 1:100, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.