



# PORCN Antibody, FITC conjugated

|                            |                                                                                                      |
|----------------------------|------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA887958LC01HU                                                                                   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                        |
| <b>Uniprot No.</b>         | Q9H237                                                                                               |
| <b>Immunogen</b>           | Recombinant Human Protein-serine O-palmitoleoyltransferase porcupine protein (269-332AA)             |
| <b>Raised In</b>           | Rabbit                                                                                               |
| <b>Species Reactivity</b>  | Human                                                                                                |
| <b>Form</b>                | Liquid                                                                                               |
| <b>Conjugate</b>           | FITC                                                                                                 |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4                     |
| <b>Purification Method</b> | >95%, Protein G purified                                                                             |
| <b>Isotype</b>             | IgG                                                                                                  |
| <b>Clonality</b>           | Polyclonal                                                                                           |
| <b>Alias</b>               | Protein-serine O-palmitoleoyltransferase porcupine (EC 2.3.1.-) (Protein MG61), PORCN, MG61 PORC PPN |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                 |
| <b>Research Area</b>       | Tags & Cell Markers                                                                                  |
| <b>Target Names</b>        | PORCN                                                                                                |