



PLPP2 Antibody, FITC conjugated

Product Code	CSB-PA018416LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43688
Immunogen	Recombinant Human Phospholipid phosphatase 2 protein (109-162AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phospholipid phosphatase 2 (EC 3.1.3.4) (Lipid phosphate phosphohydrolase 2) (PAP2-gamma) (PAP2-G) (Phosphatidate phosphohydrolase type 2c) (Phosphatidic acid phosphatase 2c) (PAP-2c) (PAP2c), PLPP2, LPP2 PPAP2C
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PLPP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.